

# KPL Milk Diluent/Blocking Solution Concentrate

Catalog No.	<u>Size</u>
5140-0011 (50-82-01)	200 mL
5140-0012 (50-82-02)	1000 mL

#### **DESCRIPTION**

KPL Milk Diluent/Blocking Solution Concentrate serves a dual purpose. As a blocking solution it binds to any remaining protein binding sites after the primary antigen or antibody is bound to the solid phase surface. As a diluent for antigens, antibodies, and other proteins, it ensures minimal non-specific interaction of these reactants with each other and with the solid phase surface.

#### FORM/STORAGE/STABILITY

5140-0011 (50-82-01) contains 2 x 100 mL 5140-0012 (50-82-02) contains 1 x 1000 mL

Store refrigerated at 2-8 °C. Stable for a minimum of one year from date of receipt when stored at 2-8 °C.

#### CONTENT

KPL Milk Diluent/Blocking Solution Concentrate contains 2% nonfat dry milk in a borate buffer.

# USE

## **Preparation:**

Before use, dilute with reagent quality water as follows. Prepare enough blocking solution to completely cover the treated surface.

	Blocking	Diluent
ELISA	1:20 (i.e. 1 mL Milk Diluent/Block + 19 mL water)	1:20 (i.e. 1 mL Milk Diluent/Block + 19 mL water)
Blotting	1:10 (i.e. 1 mL Milk Diluent/Block + 9 mL water)	1:20 (i.e. 1 mL Milk Diluent/Block + 19 mL water)

## Application:

Use the blocking solution to treat membrane or plastic surfaces **after** they have been sensitized with primary antibody or antigen. Use as a diluent for antigens, antibodies and other proteins only after sensitization of the microwell plate or membrane.

#### **Incubation Time:**

Block for at least 10 minutes to ensure adequate blocking. Proceed with assay following standard procedures.

## PRODUCT SAFETY AND HANDLING

See MSDS (Material Safety Data Sheet) for this product.

RELATED PRODUCTS	Cat. No.
Coating Solution Concentrate	5150-0014
	(50-84-00)
Wash Solution Concentrate	5150-0008
	(50-63-00)
BSA Diluent/Blocking Solution	5140-0006 (50-61-00)

508.244.6400 • 800.676.1881 Toll Free • 508.634.3334 Fax www.seracare.com